GRADIA® DIRECT Flo

LIGHT-CURED FLOWABLE COMPOSITE

GRADIA® DIRECT Flo is a multi-purpose light-cured, radiopaque fluoride releasing, micro-filled hybrid type composite resin with excellent workflow. GRADIA® DIRECT Flo allows for effective direct injection using a special dispensing tip, resulting in aesthetic restorations.

For use only by a dental professional in the recommended indications.

RECOMMENDED INDICATIONS

1. Restoration of Class I, II, III, IV, V cavities (particularly for small Class I cavities / shallow Class V cavities / other small cavities).
2. Restoration of root surface caries.
3. Restorations in deciduous teeth.
4. Filling tunnel shaped cavities.
5. Sealing hypersensitive areas.
6. Liner / base / filling in cavity undercuts.
7. Sealant.
8. Splinting mobile teeth.
9. Additions to composite restorations. (Fig. 1)

CONTRAINDICATIONS

1. Pulp capping.
2. In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

DIRECTIONS FOR USE

1. Preparations
Hold the syringe upright and remove the wing cap by turning counter-clockwise. Take care not to expose material to direct light from the dental lamp or natural light (Fig. 2).
2. Promptly and securely attach the dispensing tip (plastic or needle type) to the syringe by turning clockwise (Fig. 2)(Fig. 3).

Note :
Take care not to attach the dispensing tip too tightly. This may damage it's screw.

3. After attaching the dispensing tip, protect it with the cover until ready to use in order to avoid exposure to light (Fig. 2).

4. Cav Preparation
Prepare cavity using standard techniques.

• For pulp capping use calcium hydroxide.

5. Shade Selection
Select shade from 7 shades of A1, A2, A3, A5, A3, CV and BW.

A shades are based on Vita® Shade.

*Vita® is a registered trademark of Vita Zahnfabrik, Bad Säckingen, Germany.

4. Use of Light-cured Bonding System

For bonding GRADIA® DIRECT Flo to tooth structure, use a light-cured bonding system (Fig. 4).

G-ænial® Bond, G-Premio BOND or G-BOND which contains adhesive monomer (4-MET) for superior bondability to the structure is recommended.

Note :

When using light-cured bonding system, follow manufacturer's instructions for use.

5. Placement of GRADIA® DIRECT Flo

1) Remove cover from the dispensing tip on the syringe (Fig. 5).

2) Prior to extrusion of material, gently trial push the syringe plunger outside the mouth to make sure that the dispensing tip is securely attached to the syringe. To remove any air from the dispensing tip, with the tip pointing upwards gently push forward the syringe plunger until material reaches the mouth of the tip (Fig. 6).

Note :
If there is air inside the dispensing tip, air bubbles may be formed at the time of injection.

3) Place the dispensing tip as close as possible to the cavity, and slowly push the plunger to inject material into it (Fig. 7). Alternatively, dispense material onto a mixing pad and transfer to the cavity using a suitable instrument.

Note :
a. When attaching the dispensing tip, make sure that no material is sticking to the joint between the tip and the syringe in order to ensure a tight connection.

b. If the syringe does not extrude smoothly, remove the dispensing tip and extrude material directly from the syringe to make sure that material is coming out.

c. The material will start to harden if exposed to the dental light or ambient light. Be sure to protect it from light when working for a mixing pad.

d. After use, immediately remove the dispensing tip and tightly close the syringe with the wing cap.

4. Cleaning and Disinfecting:
MULTI-USE DELIVERY SYSTEMS: to avoid cross-contamination between patients this device requires mid-level disinfection. Immediately after use inspect device and for deterioration. Discard device if damaged.

DO NOT IMMERSO. Thoroughly clean device to prevent drying and accumulation of contaminants. Disinfect with a mid-level registered healthcare-grade infection control product according to regional/national guidelines.

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Clinical Hint 1

In order to inject effectively, use the surface tension of the material to ensure uniformity across the entire surface of the restoration during build up. Once the required amount has been injected, release the pressure on the plunger and withdraw the syringe in a direction perpendicular to the surface. This will allow the material to separate from the dispensing tip and provides a smooth surface over the restoration.

Clinical Hint 2

When filling a large cavity, it is recommended to place material incrementally into the cavity. Another effective method is to use GRADIA® DIRECT Flo for filling in undercuts or as a liner/base, and then to place composite resin (GRADIA® DIRECT, G-ænial®, Essentia®, etc.) on top.

6. Light Curing
Light the GRADIA® DIRECT Flo using a light curing unit (Fig. 8). Refer to the following chart for Irradiation Time and Effective Depth of Cure.

Note :
When light curing material, wear protective glasses.

Irradiation Time and Effective Depth of Cure

	Shade	A1, A2, A3	A3.5, BW	A03, CV
20 sec. Halogen / LED (700mW/cm²)	2.0 mm	1.5 mm	1.5 mm	
10 sec. High power LED (more than 1200mW/cm², e.g. GC-D-Light Pro)				
40 sec. Halogen / LED (700mW/cm²)	3.0 mm	2.5 mm	3.0 mm	

7. Shaping and Polishing

Shape and polish using standard techniques.

SHADES

A1, A2, A3, A5, A03, CV (Cervical color), BW (Bleaching White)

STORAGE

Recommended for optimal performance, store in a cool and dark place (4 - 25°C / 39.2 - 77.0°F)

(shelf life : 3 years from date of manufacture)

PACKAGES

1. Syringe 1.5g (0.8ml) x 2, 4 dispensing tips (needle type), 2 light protective covers
2. Dispensing tip package:
 - 20 dispensing tips (needle type), 2 light protective covers
 - 20 dispensing tips (plastic type), 2 light protective covers

CAUTION

1. In case of contact with oral tissue or skin, remove immediately with cotton or a sponge soaked in alcohol. Flush with water.

2. In case of contact with eyes, flush immediately with water and seek medical attention.

3. Take special care that the patient should not swallow material.

4. Do not mix with other similar products.

5. The dispensing tip cannot be sterilized in an autoclave or chemilavage.

6. Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.

Some products referenced in the present IFU may be classified as hazardous according to GHS.

Always familiarize yourself with the safety data sheets available at:

<http://www.geurope.com>

They can also be obtained from your supplier.

CLEANING AND DISINFECTION:

MULTI-USE DELIVERY SYSTEMS: to avoid cross-contamination between patients this device requires mid-level disinfection. Immediately after use inspect device and for deterioration. Discard device if damaged.

DO NOT IMMERSO. Thoroughly clean device to prevent drying and accumulation of contaminants.

Disinfect with a mid-level registered healthcare-grade infection control product according to regional/national guidelines.

Note :
If there is air inside the dispensing tip, air bubbles may be formed at the time of injection.

3) Place the dispensing tip as close as possible to the cavity, and slowly push the plunger to inject material into it (Fig. 7). Alternatively, dispense material onto a mixing pad and transfer to the cavity using a suitable instrument.

Note :
a. When attaching the dispensing tip, make sure that no material is sticking to the joint between the tip and the syringe in order to ensure a tight connection.

b. If the syringe does not extrude smoothly, remove the dispensing tip and extrude material directly from the syringe to make sure that material is coming out.

c. The material will start to harden if exposed to the dental light or ambient light. Be sure to protect it from light when working for a mixing pad.

d. After use, immediately remove the dispensing tip and tightly close the syringe with the wing cap.

4. Cleaning and Disinfecting:
MULTI-USE DELIVERY SYSTEMS: to avoid cross-contamination between patients this device requires mid-level disinfection. Immediately after use inspect device and for deterioration. Discard device if damaged.

DO NOT IMMERSO. Thoroughly clean device to prevent drying and accumulation of contaminants.

Disinfect with a mid-level registered healthcare-grade infection control product according to regional/national guidelines.

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