Prior to use, please read the instructions for use. EN

OPTIGLAZE color

LIGHT-CURED CHARACTERIZATION MATERIAL.

For use only by a dental professional in the recommended indications.

RECOMMENDED INDICATIONS

- Characterization of direct & indirect composite restorations, acrylic denture base and artificial acrylic teeth.
- For obtaining surface smoothness and wear resistance of restorations made of composite resin, acrylic denture base and artificial acrylic teeth.

Note:
Not intended for characterization of natural teeth.

CONTRAINDICATIONS

In rare cases the product may cause sensitivity to some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

DIRECTIONS FOR USE

Surface characterization
1. After contouring the resin surface, finish with a carbide bur or coarse silicone point. The oxygen inhibition layer must be removed as it inhibits the polymerization of OPTIGLAZE color.
Laboratory use:
Saturate with 25-50 µm alumina (0.15 MPa/1.5 bar), clean and dry.
Note:
1) If surface still contains oil residue, the surface can be cleaned further with ethanol (alcohol). The use of ethanol (alcohol) on acrylic resin may cause micro cracks. A mild detergent can be used to clean the acrylic resin.
2) When applying OPTIGLAZE color on hybrid ceramic blocks such as CERASMART, treat the surface using CERAMIC PRIMER II.

- Shake the bottle of OPTIGLAZE color well. Dispense a few drops in a dispensing dish. Dilute by adding CLEAR to the desired shade. Apply a thin layer to the resin surface using the brush provided. Do not air blow.
 - Light cure with a suitable light curing device. (See table 1.) Place the light tip as close as possible to the surface.
- Surface use:
Set a 2 cm platform in the curing device and place the resin on it. Apply the coated surface facing the light. Reverse the prosthesis and repeat the light curing to fully cure the undercut or area in shadow.

Internal characterization

- Build up the dentine layer with composite and light cure. Follow the instructions of your composite system.
Note:
Do not remove the oxygen inhibition layer. Apply OPTIGLAZE color on the oxygen inhibition layer. If the oxygen inhibition layer is removed, treat the surface with CERAMIC PRIMER II (blow dry) and COMPOSITE PRIMER I (light cure) before applying OPTIGLAZE color.
- Shake the bottle of OPTIGLAZE color well. Dispense a few drops in a disposable dispensing dish. Apply a thin layer of OPTIGLAZE color of choice to the dentine surface with an instrument or brush (e.g. the convex and concave of dentin mamelon or incisal area). Mix multiple colors to obtain the desired shade if needed. Do not air blow.
- Light cure with a suitable light curing device. (See table 1)
- Continue building up with enamel shade on the oxygen inhibition layer. Follow the instructions of your composite system.
- Finish and polish using standard technique or apply OPTIGLAZE color CLEAR or CLEAR HV on the entire surface.

Repair of restoration (e.g. when the gloss is lost)
1) Roughen the surface of restoration using a carbide bur and paper cone.
2) Clean the surface with a steam cleaner or ultrasonic cleaner and then dry with oil free air.

Note:

If surface is contaminated by oil residue, the surface can be cleaned further with ethanol (alcohol). The use of ethanol (alcohol) on acrylic resin may cause micro cracks. A mild detergent can be used to clean the acrylic resin.
1) Apply OPTIGLAZE color CLEAR or CLEAR HV shade to the surface. Do not air blow.
4) Light cure with a suitable light curing device. (See table 1)

Table 1: Light Curing Time	
GC LABOUGHT DUO	90 sec.
GC LABOUGHT LV-III_II	5 min.
GC STEPLIGHT SL-I	10 sec.
	Distance from light source <3 cm
	Distance from light source >3 cm
Halogen Light	20 sec.
Plasma Arc	40 sec.
LED Light (with light source having wavelength lower than 430 nm)	8 sec.
	40 sec.
Table 2: Depth of Cure	
Orange, Yellow, Redbrown, Red	0.07 mm
A Plus, B Plus, C Plus, White, Ivory White, Pink Orange, Pink, Olive, Lavender, Grey, Blue	0.15 mm
Clear, Clear HV	0.5 mm

STORAGE
Recommended for optimal performance, store at room temperature of (4-25°C) (39-77.0°F).
Shelf life: 3 years from date of manufacture

PACKAGE:

Complete set:
A Plus 2.6mL, B Plus 2.6mL, C Plus 2.6mL, White 2.6mL, Ivory White 2.6mL, Yellow 2.6mL, Orange 2.6mL, Pink Orange 2.6mL, Pink 2.6mL, Red Brown 2.6mL, Olive 2.6mL, Lavender 2.6mL, Grey 2.6mL, Blue 2.6mL, Red 2.6mL, Clear 5mL, Clear HV 5mL, Disposable Dispensing Dish, Flat Brush, Round Brush

CAUTION

- Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.
- In particular avoid contact of uncured material with skin.
- Avoid getting material on clothing.
- OPTIGLAZE color is highly flammable liquid and vapor. Do not use near naked flame. Keep away from source of ignition.
- OPTIGLAZE color is volatile and may cause respiratory irritation. Use in a well ventilated space. Replace cap immediately after dispensing.
- It is not recommended to use OPTIGLAZE color on high temperature burning area.
- In case of contact with eyes, flush immediately with water and seek medical attention.
- In case of contact with skin, immediately remove with a sponge or cotton pellet. Flush with water.
- The oxygen inhibited layer of the resin surface to be coated shall be removed before applying OPTIGLAZE color as it may hinder the proper curing of OPTIGLAZE color.
- Dispense a few drops of OPTIGLAZE color into a disposable dispensing dish immediately before applying as the material contains a volatile solvent.
- Store locked up.
- Dispose all waste according to local regulation.

CLEANING AND DISINFECTING
MULTI-USE DELIVERY SYSTEMS: to avoid cross-contamination
between patients this device requires mid-level disinfection. Immediately after use inspect device and label for deterioration. Discard device if damaged.
DO NOT IMMERGE. Thoroughly clean device to prevent drying and accumulation of contaminants. Disinfect with a mid-level registered healthcare-grade infection control product according to regional/national guidelines.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at:
http://www.gceurope.com or for The Americas
http://www.gcamerica.com
They can also be obtained from your supplier.

Last revised 06/2016

Rx Only

CE 0086

IGCI

MANUFACTURED by
GC DENTAL PRODUCTS CORP.
2-285 Tonimatsu-cho, Kasugai, Aichi 486-0844, Japan
RESPONSIBLE MANUFACTURER IN CANADA:
GC AMERICA INC.
3717 West 127th Street, Alsip, IL 60803 U.S.A.
TEL: +1 708 597 0000

DISTRIBUTED BY
GC CORPORATION
76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan
GC EUROPE N.V.
Reserchpark Hasroede-Leuven 1240, Interfeuvlaan 33,
B-3001 Leuven, Belgium TEL: +32 16 74 10 10
GC AMERICA INC.
3717 West 127th Street, Alsip, IL 60803 U.S.A.
TEL: +1 708 597 0000

GC ASIA DENTAL PTE. LTD.
11 Tampines Concourse, #03-05, Singapore 528729
TEL: +65 6546 7588
GC AUSTRALASIA DENTAL PTY. LTD.
17/33 Balfour Rd, Balmeadow, NSW 2019, Australia
TEL: +61 2 9301 8200

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