Saliva-Check BUFFER

In Vitro Test for checking the Quality, pH and Buffering Capacity of Saliva

For use only by a dental professional in the recommended indications.

RECOMMENDED INDICATIONS

To check the properties of resting and stimulated saliva.

CONTRAINDICATIONS

In vivo diagnostic use.

DIRECTIONS FOR USE

Prior to any visit where a saliva diagnostic procedure is planned, instruct the patient not to smoke, consume food or drink, brush the teeth or use a mouth wash for at least one hour prior to the scheduled appointment time.

I. TESTING OF RESTING SALIVA

TEST 1 – Visual inspection of level of hydration Visually assess the lower lip labial gland secretion. Evert the lower lip, gently blot the labial mucosa with a small piece of gauze and observe the mucosa under good light. Droplets of saliva will form at the orifices of the minor glands.

Assess the time for visible production of saliva as follows¹:

Greater than 60 seconds : resting flow Low ●
Less than 60 seconds : resting flow Normal ●

TEST 2 – Saliva consistency

Visually assess the resting salivary consistency in the oral cavity¹.

Sticky frothy saliva residues: Increased viscosity ●
Frothy bubbly saliva: Increased viscosity ●
Watery clear saliva: Normal viscosity ●

TEST 3 – pH measurement

Instruct the patient to expectorate any pooled saliva into the collection cup. Take a pH test strip, place this into the sample of resting saliva for 10 seconds, and then check the colour of the strip. This should be compared with the testing chart available in the package.

II. TESTING OF STIMULATED SALIVA

TEST 4 – Saliva quantity

Instruct the patient to chew the piece of wax to stimulate salivary flow. After 30 seconds, let the patient expectorate into the spittoon. Continue chewing for a further 5 minutes, collecting all the saliva into the collection cup at regular intervals.

The quantity of saliva can be measured by checking the mL markings on the side of the cup.

Quantity of saliva at 5 minutes²

< 3.5 mL Very low ●
Between 5.0 – 3.5 mL Low ●
> 5.0 mL Normal ●

Note:

Normal stimulated saliva flow rate may vary between 1mL/min – 1.6mL/min.

TEST 5 – Buffering capacity


- Remove a Buffer test strip from the foil package and place onto an absorbent tissue with the test side up.
- Using a pipette, draw sufficient saliva from the collection cup and dispense one drop onto each of the 3 test pads. Immediately turn the strip 90° to soak up excess saliva on the absorbent tissue. This will prevent the excess saliva from swelling on the test pad and possibly affecting the accuracy of the test result.
- The test pads will begin to change colour immediately and after 2 minutes the final result can be calculated by adding the points according to the final colour of each pad. See conversion table and examples underneath.

Conversion table

Test pad colour at 2 minutes

Green	4 points
Green/Blue	3 points*
Blue	2 points
Red/Blue	1 point *
Red	0 points

Examples

Interpreting the result.

Combined total	Buffering ability of Saliva
0-5	Very low ●
6-9	Low ●
10-12	Normal / High ●

STORAGE

Recommended for optimal performance, store components at room temperatures (15°-30°C). (Shelf life: 3 years from date of manufacture)

PACKAGE

Contents	Quantity
In vitro pH test strips	20
Saliva collection cup	20
Wax for saliva stimulation	20
Saliva dispensing pipette	20
Buffer test strip	20

CAUTION

- In warm ambient temperatures, store the wax in a refrigerator.
- After opening the wax and Buffer test strip, use immediately.
- All components used for each test are for single use only. Do not reuse.
- Dispose of according to local regulations.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at: <http://www.gceurope.com> or for The Americas <http://www.gcamerica.com>. They can also be obtained from your supplier.

¹ Information by kind permission of Prof. L Walsh©2002

² Lund University Faculty of Odontology Department of Cariology, Sweden

Last revised: 01/2018

Rx Only

DISTRIBUTED by
GC CORPORATION
76-1 Hatanuma-cho, Itabashi-ku, Tokyo 174-8585, Japan

EU: GC EUROPE NV
Ruessepark Hassode-Leuven 1240, Interleuvenaan 33,
B-3001 Leuven, Belgium TEL: +32 16 74 10 00

GC AMERICA INC.
3737 West 127th Street, Alsip, IL 60803 U.S.A.
TEL: +1 708-597-0900
www.gcamerica.com

GC SOUTH AMERICA
Rua Heliodora, 399, Santana - São Paulo, SP, BRASIL
CEP: 02022-091 - TEL: +55-11-2925-0965
CNPJ: 06.278.999/0001-61

RESP: TEC: Mayara de Santa Ribeiro - CRO/SP 105.882
GC ASIA DENTAL PTE. LTD.
11 Tampines Concourse, #03-05 Singapore 528729
TEL: +65 6546 7988

GC AUSTRALASIA DENTAL PTY. LTD.
1753 Botany Rd, Banksmeadow, NSW 2019, Australia
TEL: +61 2 9301 8200

MANUFACTURED by
SHENZHEN KANG SHENG BAO
BIO-TECHNOLOGY CO., LTD.

Second Floor, Block 3, An Da Electronic Industrial Park, Fu Hai,
Bao An District, 518103 Shenzhen, Guangdong Province, PRC.

RESPONSIBLE MANUFACTURER IN CANADA:
GC AMERICA INC.
3737 West 127th Street, Alsip, IL 60803 U.S.A.
PRINTED IN HONGKONG

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