Prior to use, carefully read the instructions for use.

EN

G-aenial™ Universal Injectable

UNIVERSAL LIGHT-CURED RADIOPAQUE INJECTABLE COMPOSITE

For use only by a dental professional in the recommended indications.

INDICATIONS FOR USE

- Direct restorative for Class I, II, III, IV and V cavities
- Sealing hypersensitive areas
- Temporary restorations, temporary crown & bridge, defect margins when margins are in enamel
- Blocking out undercuts
- Liner or base
- Temporary crowns & bridges, inlays and veneers using the indirect technique in combination with GRADIA® or GRADIA PLUS components (please refer to their dedicated IFU).
- Splinting of teeth in combination with fibres such as GC everStick fibres. Follow manufacturer's instructions.

CORONARATIONS

- Direct pulp capping
- Avoid use of this product in patients with known allergies to methacrylate monomer or methacrylate polymer.

DIRECTIONS FOR USE

- Shade Selection:** From the 16 shades of XBW, BW, A1, A2, A3, A5, A4, B1, B2, CV, CVD, AO1, AO2, JE, AE. A, B, C, D shades are based on Vita® Shade.
- Cavity Preparation:** Prepare cavity using standard techniques. Dry by gently blowing with oil free air. Rubber dam is recommended to isolate the prepared tooth from contamination. Note: For pulp capping, use calcium hydroxide.

- For bonding:** G-aenial Universal Injectable to enamel and/or dentin, use a light-cured bonding system such as G-Premio BOND, G-aenial Bond or G-BOND. Follow manufacturer's instructions.

- For direct pulp capping:**
 - Dispense from a syringe.
 - Hold the syringe upright and remove the wing cap by turning counter-clockwise (Fig. 1).
 - Promptly and securely attach the dispensing tip to the syringe by turning clockwise.
 - Place the light protective cap onto the dispensing tip.
 - Remove cap from the dispensing tip.
 - Note: Take care not to attach the dispensing tip too tightly. This may damage its screw. In order to ensure a tight connection, make sure that the threads are free of remnant material.
 - Note: To prevent any material in the cavity, make sure to check the following:
 - Dispensing tip is securely attached to the syringe.
 - Any possible air from the dispensing tip is removed by gently pushing forward the syringe plunger with the thumb (Fig. 2).


Fig. 1

Fig. 2

Fig. 3

- For bonding, G-aenial Universal Injectable should be applied to enamel and/or dentin using a light-cured bonding system such as G-Premio BOND, G-aenial Bond or G-BOND. Follow manufacturer's instructions.
- For direct pulp capping, prepare the cavity using standard techniques. Dry by gently blowing with oil free air. Rubber dam is recommended to isolate the prepared tooth from contamination. Note: For pulp capping, use calcium hydroxide.
- For bonding: G-aenial Universal Injectable to enamel and/or dentin, use a light-cured bonding system such as G-Premio BOND, G-aenial Bond or G-BOND. Follow manufacturer's instructions.
- For direct pulp capping:
 - Dispense from a syringe.
 - Hold the syringe upright and remove the wing cap by turning counter-clockwise (Fig. 1).
 - Promptly and securely attach the dispensing tip to the syringe by turning clockwise.
 - Place the light protective cap onto the dispensing tip.
 - Note: Take care not to attach the dispensing tip too tightly. This may damage its screw. In order to ensure a tight connection, make sure that the threads are free of remnant material.
 - Note: To prevent any material in the cavity, make sure to check the following:
 - Dispensing tip is securely attached to the syringe.
 - Any possible air from the dispensing tip is removed by gently pushing forward the syringe plunger with the thumb (Fig. 2).

- Place the dispensing tip as close as possible to the prepared cavity, and slowly push the plunger to extrude material onto a mixing spatula or transfer it to the cavity using a suitable instrument.

- Note: If it plugged, remove the dispensing tip and extrude a small amount of material directly from the syringe.
- Minor exposure to ambient light. Ambient light can shorten the mixing time.
- After use, immediately remove and dispose of the dispensing tip and tightly close the syringe with the wing cap.

- For dispensing:** Insert the G-aenial Universal Injectable Unit into the Unipit APPLIER II or equivalent. Remove the cap and extrude material directly into the prepared cavity. Use steady pressure (Fig. 3). Maintain pressure on the applier until the dispensing tip is removed from the mouth. This prevents the Unipit from coming loose from the applier.

- Note: a. The material may be difficult to extrude immediately after removing from the applier. Prior to use, wait for a few minutes at room temperature.
- b. For first use, apply pressure slowly and gradually to control extrusion of material.

- Curing:** In order to inject effectively, use the surface tension of the material to ensure uniformity across the entire surface of the restoration during build up. Once the required amount has been injected, release the pressure on the plunger and withdraw syringe in a direction perpendicular to the surface. This will allow the dispensing tip to separate from the dispensing tip and provides a smooth surface over the restoration.

- Light-curing:** Light-cure G-aenial Universal Injectable using a light curing unit. Keep light guide as close as possible to the surface. Refer to the following chart for Irradiation Time and Effective Depth of Cure (ISO 4049)

	10 sec. (High-intensity LED) (>1200 mW/cm²)	20 sec. (Halogen/ED) (>700 mW/cm²)
A1, A2, A3, B1, B2, JE, AE	2.5 mm	2.0 mm
XBW, BW, A3, 5, A4, AO1, AO2, CV, CVD	2.0 mm	1.5 mm
AO3		

- Note: a. material should be placed and light-cured in layers. For maximum layer thickness, please consult above table.
- b. Lower light intensity may cause insufficient curing or discoloration of the material.
- c. For the indirect technique, light-cure each layer for 30 seconds with LABOLIGHT LV-II. If for 10 seconds with STEPLIGHT SL-I. Final cure for 3 minutes using LABOLIGHT LV-II.

- Finishing and Polishing:** Finish and polish using standard techniques.

- SHADES:** XBW, BW, A1, A2, A3, A5, A4, B1, B2, CV, CVD, AO1, AO2, AO3, JE, AE

- STORAGE:** Recommended for optimal performance, store in a cool and dark place (4°C-25°C/9-27°F) away from high temperatures or direct sunlight.

- PACKAGES:**
 - I. Syringes
 - 1 Syringe 1.7 g (0.16 mL) x 1, Dispensing Tip Long Needle x 10, Light-protective cap x 10
 - 2 Syringe 1.7 g (0.16 mL) x 1, Dispensing Tip Needle x 5
 - 3 Syringe 1.7 g (0.16 mL) x 2, Dispensing Tip Needle x 10
 - 4 Dispensing Tip Refill
 - 5 Dispensing Long Needle x 30, Light-protective cap x 3
 - 6 Finishing and Polishing

- II. Units

- Units 0.27 g (0.16 mL) x 15

- (All packaging types may not be available in all countries.)

- CAUTION:**
 - In case of contact with oral tissue or skin, remove immediately with cotton or a sponge soaked in alcohol. Flush with water. To avoid contact, a rubber dam and COCOA BUTTER can be used to isolate the operation field from the skin or oral cavity.
 - In case of contact with eyes, flush immediately with water and seek medical attention.
 - To prevent ingestion of the material.
 - After starting the filling process, partially extract paste outside of the patient's mouth and away from the patient to ensure paste is flowing properly.
 - Do not mix with other dental products.
 - Do not use for single use only. To prevent cross contamination between patients, do not reuse the dispensing tip. The dispensing tip cannot be sterilized in an autoclave or chemosteamer.
 - For single use, wear plastic or rubber gloves during operation to avoid direct contact with air inhibited resin layers in order to prevent single use only.
 - Wear protective eye glasses during light curing.
 - Wear protective clothing during cutting.
 - Not to mask or cover the nose with a dust collector and wear a dust mask to avoid inhalation of cutting dust.
 - Avoid getting material on clothing.
 - In case of contact with any dental or ortho appliance, spoon or cotton pellet before light-curing.
 - Do not use G-aenial Universal Injectable in combination with eugenol containing materials as eugenol may hinder G-aenial Universal Injectable.
 - Do not wipe the dispensing tip with excessive alcohol as this may cause contamination of the material.
 - Dispose of all wastes according to local regulations.
 - Some products referred to in the present IFU may be classified as hazardous according to GHS. Always familiarise yourself with the safety data sheets available <http://www.europeure.com>.

- The material should also be obtained from your supplier.

- CLEANING AND DISINFECTION:** MULTI-USE DELIVERY SYSTEMS: to avoid cross-contamination between patients this device requires mid-level disinfection. Immediately after use inspect device and label for deterioration. Discard if deteriorated.

- DO NOT USE: G-aenial Universal Injectable with a mid-level registered healthcare-grade infection control product according to regional / national guidelines.

- For the Summary of Safety and Clinical Performance (SSCP) please see EU DAMED database (<https://ec.europa.eu/tools/eudamed>) or contact us at Regulatory.gce@gc.dental

- Unpublished effects- Reporting:

- If you become aware of any kind of undesired effect, reaction or similar events experienced by use of this product, including those not listed in this instruction for use, please report them directly through the relevant vigilance system, by selecting the proper authority from your country accessible through the following link: https://ec.europa.eu/growth/sectors/medical-devices/contacts_en, as well as to our internal vigilance system: <vigilance@gc.dental>. In this way you will contribute to improve the safety of this product.

- GC UNITED KINGDOM Ltd. Coopers Court Newport Pagnell, Buckinghamshire, MK16 8JS, United Kingdom

- Last revised: 10/2022

Prior to use, carefully read the instructions for use.

EN

Prior to use, carefully read the instructions for use.

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