Prior to use, carefully read
the instructions for use.

EN

GC Fuji® II LC CAPSULE

RADIOPAQUE LIGHT CURED REINFORCED
GLASS IONOMER RESTORATIVE IN CAPSULES

For use only by a dental professional in the
recommended indications.

RECOMMENDED INDICATIONS

1. Restoration of Class III, V and limited Class I
cavities.
2. Restoration of primary teeth.
3. Core build up.
4. Cases where a radiopaque restoration is required.
5. Geriatric applications.
6. As a base or liner.

CONTRAINdications

1. Pulp capping.
2. In rare cases the product may cause sensitivity in
some people. If any such reactions are experienced,
discontinue the use of the product and refer to a
physician.

DIRECTIONS FOR USE

Powder / Liquid Ratio (g/g) 0.33 / 0.10

Mixing Time (sec.) 10"

Working Time (min., sec.) 3'15"

Light Curing Time (sec.) 20"

Depth of Cure (A2) (mm) 1.8

Test conditions : Temperature (23 +/- 1°C)
Relative humidity (50 +/- 10%)

1. CAVITY PREPARATION

- a) Prepare tooth using standard techniques.
Extensive mechanical retention is not necessary.
For pulp capping, use calcium hydroxide.
- b) Wash and dry but DO NOT DESICCATE.
GC CAVITY CONDITIONER or GC DENTIN
CONDITIONER is recommended to remove the
smear layer (Fig. 1).

2. CAPSULE ACTIVATION AND MIXING

- a) Before activation, shake the capsule or tap its
side on a hard surface to loosen the powder
(Fig. 2).
- b) To activate the capsule, push the plunger until it
is flush with main body (Fig. 3).
- c) Immediately place the capsule into a metal
GC Capsule Applier and click the lever once
(Fig. 4). The capsule is now activated.
Note :
The capsule should be activated just before
mixing and used immediately.
- d) Immediately remove the capsule and set it into a
mixer (an amalgamator) and mix for 10 seconds
at high speed (+/-4,000 RPM) (Fig. 5).

3. RESTORATIVE TECHNIQUE

- a) Immediately remove the mixed capsule from the
mixer and load it into the GC Capsule Applier.
 - b) Make two clicks to prime the capsule then
syringe (Fig. 6).
- The working time is 3 minute 15 seconds from
start of mixing at 23°C (73.4°F). Higher
temperature will shorten working time.
Note :
1) To adjust the direction of the nozzle, hold the
applier with the capsule towards you and turn
the capsule body.
2) To remove the used capsule, push the applier
release button. Twist the capsule and pull
upwards.
- c) Remove surface moisture but DO NOT
DESICCATE.
- d) Extrude cement directly into preparation (Fig. 7).
Avoid air bubbles.
- e) Form the contour and place a matrix if required.
- f) Light-cure for 20 seconds using a suitable visible
light curing device (470nm wavelength) (Fig. 8).
Place light source as closely as possible to the
cement surface.
- Note :
For cavities deeper than 1.8 mm, use a layering
technique.
4. FINISHING
- Finish under water spray using superfine diamond
bur, silicone point and polishing strips (Fig. 9).
- Note :
Apply a final coat of GC Fuji VARNISH (blow dry)
or GC Fuji COAT LC (light cure for 10 sec.) or
G-COAT PLUS (light cure for 20 sec.) (Fig. 10).

STORAGE

Recommended for optimal performance, store in a
cool and dark place (4-25°C / 39.2-77.0°F).
(Shelf life : 2 years from date of manufacture)

SHADE

A1, A2, A3, A3.5, A4, B2, B3, B4, C2, C4, D2
(11 shades based on Vita® shades).

Vita® is a registered trademark of Vita Zahnfabrik, Bad Säckingen, Germany.

PACKAGES

1. Single shade package : 50 capsules.
2. Assortment package : 50 capsules (10 pcs each of
A4, B2, B3, C2; 5 pcs each of B4, D2).
* Average contents per capsule : 0.33g powder and
0.10g (0.085mL) liquid.
- * Minimum net volume of mixed cement per capsule
: 0.10mL.

Option : GC Capsule Applier (1 piece).

CAUTION

1. In case of contact with oral tissue or skin, remove
immediately with a sponge or cotton soaked in
alcohol. Flush with water.
2. In case of contact with eyes, flush immediately with
water and seek medical attention.
3. Personal protective equipment (PPE) such as
gloves, face masks and safety eyewear should
always be worn.

Some products referenced in the present IFU may
be classified as hazardous according to GHS.
Always familiarize yourself with the safety data
sheets available at: <http://www.gceurope.com>
They can also be obtained from your supplier.

Last revised : 11/2016

CE 0086

'GC'

MANUFACTURED by
GC CORPORATION
76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan

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