

GC Fuji BOND LC

LIGHT CURED GLASS IONOMER BONDING AGENT



**WE DON'T SHARE THE OPINION
ABOUT TOTALLY SEALED
COMPOSITE FILLINGS**


Masticatory pressures weaken the marginal seal of composite fillings in the long term and help to bring about marginal leakage. In spite of extensive research, even modern adhesives do not possess adequate elasticity, or properties effective in long-term caries prevention. Marginal leakage results in secondary caries.

The search for a material which could cushion masticatory pressures - like a shock-absorber - and ensure adequate protection against caries led to the development of **Fuji BOND LC**, the first modern glass ionomer - based, dentine - enamel adhesive.

Fuji BOND LC is ideally suited for surgery use, both in meeting today' s demanding clinical requirements and in its ease of use.

FUJI BOND LC, THE UNIVERSAL ADHESIVE FOR CARIES-FREE COMPOSITE FILLINGS

SIMPLE APPLICATION



Its ease of application is impressive. For the first time, a glass ionomer cement of free-flowing consistency can be very rapidly mixed with a brush.




Apply **GC Cavity Conditioner** to the cavity surfaces, leave for 10 sec., rinse off and dry. **Cavity Conditioner** removes the smearlayer and seals the dentinal tubules, thereby greatly enhancing adhesion and protection against sensitivity.



Immediately apply **Fuji BOND LC** with a brush, like existing adhesives, and light-cure for 20 seconds. Apply composite, light-cure. Ready. Preparation of the cement borders is unnecessary.



SAFETY PACKAGE



The safety package to prevent secondary caries is equally impressive. Prolonged fluoride release into the surrounding tooth structure provides effective protection against secondary decay.

■ Fuji BOND LC - adhesive
■ Fuji B LC - restorative
■ Fuji LUTE - luting cement

Compensation for shrinkage




Compensation for masticat. pressure



Thanks to its high elasticity, **Fuji BOND LC** continually cushions masticatory pressures (like a shock-absorber), without disturbing the bond between the tooth structure and the glass ionomer filling.

In summary, **Fuji BOND LC** offers a universal adhesive combining impressive, practice-oriented manipulation times with the well-known advantages of glass ionomer materials.



PACKAGING INFORMATION

- STANDARD SET
 - 1 x 5 g powder
 - 1 x 7 ml liquid
 - 1 x 5.7 ml Cavity Conditioner
- REFILLS
 - 5 g powder
 - 7 ml liquid
 - 5.7 ml Cavity Conditioner

GC EUROPE N.V. - Head Office
 Interleuvenlaan 13, B - 3001 Leuven,
 Tel. +32.16.39.80.50, Fax. +32.16.40.02.14
 E-mail info@gceurope.com, http://www.gceurope.com

GC UNITED KINGDOM Ltd.
 22, Coopers Court, Newport Pagnell
 UK - Bucks. MK16 8JS
 Tel. +44.1908.218.999, Fax. +44.1908.218.900
 E-mail info@uk.gceurope.com, http://uk.gceurope.com