GC Gold Label Light-Cured Universal Restorative

RADIOPAQUE LIGHT CURED REINFORCED GLASS IONOMER RESTORATIVE

For use only by a dental professional in the recommended indications.

RECOMMENDED INDICATIONS

- Class III and V restorations, particularly for cervical erosions and root surface caries.
- Restoration of primary teeth.
- Core build up.
- Cases where radiopacity is required.
- Geriatric applications.
- As a base or liner.

CONTRAINDICATIONS

- Pulp capping.
- In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

DIRECTIONS FOR USE

Powder / Liquid Ratio (g/g)	3,2 / 1,0
Mixing Time (sec.)	20"-25"
Working Time (min., sec.)	3'45"
Light Curing Time (sec.)	20"
Depth of Cure (A2) (mm)	1,8

Test conditions: Temperature (23+/- 1°C) Relative humidity (50+/-10%)

- POWDER AND LIQUID DISPENSING**
 - Select shade based on Vita® shade guide.
® : Vita is a registered trademark of Vita Zahnfabrik, Bad Säckingen, Germany.
 - The standard powder to liquid ratio is 3.2g / 1.0g. 1 level scoop of powder to 2 drops of liquid.
 - For accurate dispensing of powder, tap the bottle gently. Do not shake or invert.
 - Hold the liquid bottle vertically and squeeze gently.
 - Close bottles immediately after use.
- MIXING**
 - Fluff powder of bottle before dispensing (Fig. 1).
 - Place one scoop of powder and two drops of liquid on mixing pad. Divide powder in half (Fig. 2).
 - Spread liquid out into a thin layer (about the size of a half dollar or 3 cm) with plastic spatula (Fig. 3).
 - Pull half of the powder onto liquid and mix with lapping strokes (like impression material) for 10 to 15 seconds (Fig. 4).
 - Pull in remaining powder and mix thoroughly to a glossy consistency. Do not exceed 20-25 second total mixing time (Fig. 5).
- RESTORATIVE TECHNIQUE**
 - Prepare tooth using standard techniques. Extensive mechanical retention is not necessary. For pulp capping use calcium hydroxide.
 - Apply CAVITY CONDITIONER (GC) (10 seconds) or DENTIN CONDITIONER (GC) (20 seconds) to the bonding surfaces using a cotton pellet or sponge.
 - Rinse thoroughly with water. Dry by blotting with a cotton pellet or gently blowing with an air syringe. DO NOT DESICCATE. Best results are obtained when prepared surfaces appear moist (glistening).
 - Mix the required amount of cement. At 23°C (73.4°F), working time is 3 minutes 45 seconds from the start of mixing. Higher temperatures will shorten working time.
 - Transfer cement to the prepared tooth surface using a syringe or suitable placement instrument. Avoid air bubbles.
 - Form the contour and place a transparent matrix if required.
 - Light cure for 20 seconds using a visible light curing device (470 nm wavelength). Place light source as close as possible to the cement surface.
Note :
For cavities deeper than 1.8 mm, use a layering technique.
- FINISHING**

Remove matrix and finish under water spray using standard techniques.

Note :
If a matrix has not been used, and the surface is exposed to air for more than a few minutes it should be sealed with GC Fuji VARNISH (blow dry) or GC Fuji COAT LC (light cure).

STORAGE

Recommended for optimal performance, store in a cool and dark place (4-25°C) (39.2-77.0°F).
(Shelf life : 3 years from date of manufacture.)

SHADE

(11 shades) A1, A2, A3, A3.5, A4, B2, B3, B4, C2, C4, D2.

PACKAGES

- Bottle of 15g powder with scoop (0.21g)
 - Bottle of 8g (6.8mL) liquid
 - 1-1 Package*: 15g powder, 8g (6.8mL) liquid, powder scoop (0.21g), mixing pad (No.22)
 - Mini pack*: 5g powder, 3g (2.6mL) liquid, powder scoop (0.21g)
- * 1-1 package and mini pack are only available in limited shades (1-1 package: A2, A3) (Mini pack: A2, A3, A3.5).

CAUTION

- In case of contact with oral tissue or skin, remove immediately with a sponge or cotton soaked in alcohol. Flush with water.
- In case of contact with eyes, flush immediately with water and seek medical attention.
- Do not mix powder or liquid with other glass ionomer components.
- Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.


Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available for the Americas at:
<http://www.gcamerica.com>
They can also be obtained from your supplier.

Last revised: 02/2016


MIXING PENGADUKAN
MEZCLA CAUCH TROÁN VAÄT LIEÄU
MISTURA การผสม


(Fig.1)


(Fig.2)

(Fig.3)

(Fig.4)

(Fig.5)

GC Gold Label Light-Cured Universal Restorative

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3737 West 127th Street, Alsip, IL 60803 U.S.A.
GC ASIA DENTAL PTE. LTD.