

CONTACT:
GC America Inc.
3737 W. 127th Street
Alsip, IL 60803
(708) 597-0900

FOR IMMEDIATE RELEASE:
February 7, 2005

GC GRADIA® DIRECT LoFlo with HDR & GRADIA® DIRECT Flo Flowable Composite Restoratives

GC America will be introducing at the Chicago Midwinter Meeting, GRADIA® DIRECT Flo and GRADIA® DIRECT LoFlo, light-cured flowable composites. These multi-purpose light-cured composites are both radiopaque, fluoride releasing, highly filled micro-hybrid resins. Two self-adapting formulas have been designed to best accommodate every application (technique) requirement.

GRADIA® DIRECT LoFlo takes advantage of GC's new HDR (High Density Radiopaque Pre-polymerized Filler) technology. This technology allows for excellent radiopacity while providing physical properties (strength and wear) similar to traditional composites making it ideal as a final restorative in all primary and conservative permanent restorations. Its highly filled formula ensures exceptionally low shrinkage, this along with its optimal flow characteristics, ensures better adaptation to proximal and cervical walls and internal line angels reducing the risk of contraction stress and microleakage.

GRADIA® DIRECT Flo, our less viscous version, because of its fluid handling properties is typically used to restore undercuts, as a liner/base under direct restoratives, and as a 60 second sealant when used in conjunction with GC G-BOND™ one component Self-Etching Light-Cured Adhesive.

Both products are packaged in 1.5g syringes, either 2 or 4 per box. Each box will come with a supply of 4 tips. GRADIA® DIRECT Flo and GRADIA® DIRECT LoFlo both come in 7 shades: A1, A2, A3, A3.5, AO3, CV (Cervical Color), and BW (Bleaching White).

GRADIA® DIRECT Flo and GRADIA® DIRECT LoFlo have indication for the following:

1. Restoration of Class I, II, III, IV, V cavities (particularly for small Class I cavities / shallow Class V cavities / other small cavities).
2. Restoration of root surface caries.
3. Restorations in deciduous teeth.
4. Filling tunnel shaped cavities.
5. Sealing hypersensitive areas.
6. Liner/Base/Filling in cavity undercuts.
7. Sealant.
8. Fixation of mobile teeth.
9. Additions to composite restorations.

GC Corporation is the company to trust for exciting new products with unequaled levels of quality and a total commitment to service. GC has fully dedicated research and development team committed to bringing new and innovative materials to the dental profession.

If you would like more information on GC America and its products, please visit www.gcamerica.com or call Customer Service at 800-323-7063.