



GC AMERICA INC.
 3737 W 127th STREET
 ALSIP, ILLINOIS 60803
 TEL (800) 323-3386
 FAX (708) 897-4062
 www.gcamerica.com

Contact:

Dan Boskocevic
 Dental Laboratory Product Manager
 GC America Inc.
 800-323-7063 ext. 3012

FOR IMMEDIATE RELEASE:

Alsip, IL – September 2014 - GC America Inc. introduces **OPTIGLAZE™ Color**, a light-cured characterization coating for composite and acrylic indirect restorations. **OPTIGLAZE Color** is now available for purchase.

OPTIGLAZE Color is a revolutionary new material developed specifically to stain and glaze your composite, acrylic, and PMMA type materials. You can now change shades, shift value, characterize, create effects and glaze all your resin type prosthetics in no time at all. Best of all, offer your dentists and patients a more beautiful and esthetic restoration that will last due to GC's innovative new homogeneously dispersed nano-filler technology.

The versatility of **OPTIGLAZE Color** will impress you because it adds color, surface gloss and high wear resistance to all your:

- Composite crowns and bridge
- Temporary resin crowns and bridges
- Hybrid ceramic crowns
- Implant retained dentures
- Acrylic permanently fixed with implants
- Artificial teeth
- Acrylic dentures

Characteristics and benefits:

- Ready-to-use, easy-to-handle and cuts on the polishing stage, saving you valuable time
- Choose from a wide variety of colors, both for internal and external characterization leading to perfect aesthetics
- Renowned nano-filler technology gives your restoration a high-wear resistance for long-lasting gloss with a high discoloration resistance

SKU	SKU Description	Product Description
008408	OPTIGLAZE Color	Complete set includes all colors and disposable palette, flat and round brush and brush holder
008409	OPTIGLAZE Color	A plus 2.6mL
008410	OPTIGLAZE Color	B plus 2.6mL
008411	OPTIGLAZE Color	C plus 2.6mL
008412	OPTIGLAZE Color	White 2.6mL
008413	OPTIGLAZE Color	Ivory white 2.6mL
008414	OPTIGLAZE Color	Yellow 2.6mL
008415	OPTIGLAZE Color	Orange 2.6mL
008416	OPTIGLAZE Color	Pink orange 2.6mL
008417	OPTIGLAZE Color	Red brown 2.6mL
008418	OPTIGLAZE Color	Olive 2.6mL
008419	OPTIGLAZE Color	Grey 2.6mL
008420	OPTIGLAZE Color	Blue 2.6mL
008421	OPTIGLAZE Color	Lavender 2.6mL
008422	OPTIGLAZE Color	Pink 2.6mL
008423	OPTIGLAZE Color	Red 2.6mL
008424	OPTIGLAZE Color	Clear 5mL
008425	OPTIGLAZE Color	Clear HV 5mL

###

