

Test results



Name of patient

File reference

Date

Saliva-Check Buffer / Saliva-Test

Resting saliva			Stimulated saliva		
Step 1 Hydration	Step 2 Viscosity	Step 3 pH	Step 4 Quantity	Step 5 Buffering	
> 60 sec. <input type="checkbox"/>	sticky frothy <input type="checkbox"/>	5.0-5.8 <input type="checkbox"/>	< 3.5 ml <input type="checkbox"/>	0-5 points <input type="checkbox"/>	
< 60 sec. <input type="checkbox"/>	frothy bubbly <input type="checkbox"/>	6.0-6.6 <input type="checkbox"/>	3.5-5.0 ml <input type="checkbox"/>	6-9 points <input type="checkbox"/>	
	waterly clear <input type="checkbox"/>	6.8-7.8 <input type="checkbox"/>	> 5.0 ml <input type="checkbox"/>	10-12 points <input type="checkbox"/>	

Saliva-Check Mutans / Streptococcus mutans-Test

> 500.000 cfu/ml <input type="checkbox"/>
< 500.000 cfu/ml <input type="checkbox"/>

Plaque Indicator Kit / Plaque-Test

pH	Mature plaque	Fresh plaque
5.0-5.8 <input type="checkbox"/>	visible <input type="checkbox"/>	visible <input type="checkbox"/>
6.0-6.6 <input type="checkbox"/>		
6.8-7.8 <input type="checkbox"/>	not visible <input type="checkbox"/>	not visible <input type="checkbox"/>